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
Step 1: Registration of the Sample in the Elements 4.0 Software

The sample is registered in the ELEMENTS 4.0 software. The weight is automatically transferred from the balance (step 1) – fast, safe, and error-free.

Step 2: Weighing and Application

Whether a few milligrams or up to 1 gram: the ELEMENTRAC CN-r offers the highest precision with various sample types – solid or liquid. These are simply weighed into tin foil or tin capsules and reliably analyzed.

Paso 3: Análisis

The samples are conveniently placed into the carousel. The CN-r automatically detects the carousel and clearly assigns all samples – ensuring a smooth, safe, and automated process.

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
CONFIGURACIONES

The ELEMENTRAC CN-r can be flexibly configured:


- | Nitrogen/Protein determination
- | Nitrogen/Protein and Carbon determination

While carbon is detected as CO₂, nitrogen determination is performed in elemental form using a highly sensitive thermal conductivity cell (TCD). Helium is used as the standard carrier gas (argon optional). This ensures maximum sensitivity in nitrogen analysis.

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All chemicals required for operation are ergonomically arranged at the front and safely hidden behind a door during routine use. Maintenance is especially user-friendly: combustion and reduction tubes can be replaced without tools in just a few minutes. This well-thought-out design reduces time expenditure, increases reliability, and ensures maximum reproducibility of your analyses.

TOUCHSCREEN & INTEGRATED PC

The ELEMENTRAC CN-r is equipped with a modern touchscreen that can be placed flexibly either directly on the device or conveniently on the desk. Thanks to the integrated PC, no additional hardware is required. The ELEMENTS 4.0 software can be intuitively operated via touchscreen – or classically with a keyboard and mouse. With the available USB ports, you can easily connect peripheral devices such as a keyboard, mouse, or a balance.

NETWORKING & DATA ACCESS

The CN-r can be easily integrated into your network. This allows you to transfer results directly into your LIMS or access the device software from your office PC. Whether recalculations, evaluations, or entering new samples: you always have full access and can complete your tasks easily and efficiently – regardless of the device's location.

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
- | Analysis of solid or liquid samples such as cheese, blue cheese, yogurt, and milk
- | Determination of protein content for nutritional labeling and quality control
- | The ELEMENTRAC CN-r meets or exceeds the requirements of all common international standards, such as DIN EN ISO 14891

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- | Determination of protein content for nutritional labeling and quality control
- | Verification of meat authenticity and prevention of food fraud
- | Monitoring protein levels in processed meat products

CARBÓN

- | Measurement of nitrogen content to assess quality
- | Control of nitrogen levels in coke production
- | Optimization of combustion processes and reduction of emissions

FERTILIZANTES

- | Determination of nitrogen concentration to ensure fertilizer efficiency
- | Quality control and compliance with agricultural standards
- | Evaluation of nutrient availability for improved crop yield

CEREALES & SEMILLAS OLEAGINOSAS

- | Protein determination for nutritional evaluation and feed production
- | Quality assessment in oilseed processing (soy, sunflower, rapeseed, etc.)
- | Supporting research on oilseed breeding and agricultural development

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